

Product Data Sheet

Elution Solution

Catalog Number: **ES-10, ES-100**

Revision Date: 2/19/21

For Research Use Only.

Product Description:

Elution Solution is an ala carte component of the “**Respiratory Viral RNA/DNA Extraction-Magnetic Kit**”. Elution Solution is a pH buffered solution designed to enhance release of nucleic acids (RNA and DNA) from **DNA/RNA Binding Beads (RS)** and contains chelating agents to inhibit DNASE and RNASE activity.

Contents and Storage:

Component	Catalog No.	Volume	Quantity	Storage
Elution Solution	ES-10	10 ml	1	15°C to 30°C
	ES-100	100 ml	1	

Procedure:

See “Respiratory Viral RNA/DNA Extraction-Magnetic Kit” product insert for procedural guidelines.