

Product Data Sheet

20X TBS Buffer with Tween-20

Catalog Number: **TBST-1G**

Revision Date: 2/19/21

For Research Use Only.

Product Description:

20X TBS Buffer with Tween-20 is a Tris buffered saline solution containing tween-20 (pH 7.7). The solution contains an isotonic concentration of NaCl and the surfactant Tween-20. The buffer is provided as a 20x concentrated solution and should be diluted prior to use. The diluted solution is appropriate for many immunohistochemical (IHC) procedures.

Contents and Storage:

Component	Catalog No.	Volume	Volume Reconstituted	Storage
20X TBS Buffer with Tween-20	TBST-1G	1 gallon	20 gallons	Ambient

Reagent Preparation:

- Dilute 1 part **20X TBS Buffer with Tween-20** with 19 parts distilled water.

Procedure:

- Follow your laboratories standard operating procedures where washing of histological samples is required.